

SPECIFICATION

Ethyl Acetate for Pesticide Residue And Polychlorinated Biphenyl Analysis

Concentration 300

| REQUIREMENT | SPECIFICATION |
|---|------------------------|
| * Appearance | Colorless clear liquid |
| * Suitability for pesticide residue analysis (GC-FPD) | to pass test |
| ** Assay (CH ₃ COOC ₂ H ₅) (GC) | min.99.7%(area/area) |
| Density (20°C) | 0.899~0.902g/mL |
| Water | max.0.05%(mass/mass) |
| Residue after evaporation | max.2ppm(mass/mass) |
| Suitability for pesticide residue-PCB analysis | to pass test |

JIS K 8110:2023(*:Additional test performed by Wako but not required by JIS)(**:The specification is originally set by Wako not JIS)