

Code No.05-0555

## SPECIFICATION

## 2-Ethylhexyl p-(Dimethylamino) Benzoate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ brown, clear liquid
Density (20°C)	under investigation
Refractive index 20°C	under investigation
Water	max.0.5%
Acidity [as (CH <sub>3</sub> ) <sub>2</sub> NC <sub>6</sub> H <sub>4</sub> COOH]	max.0.5%
Assay (capillary GC)	min.97.0%