

SPECIFICATION

Ethyl Acetate

for Preparative Chromatography

| REQUIREMENT | SPECIFICATION |
|-----------------------------------|------------------------|
| Appearance | Colorless clear liquid |
| Density (20°C) | 0.897~0.902g/ml |
| Refractive index 20°C | 1.372~1.374 |
| Absorbance (260nm) | max.0.40 |
| Absorbance (300~400nm) | max.0.01 |
| Water | max.0.1% |
| Residue after evaporation | max.5ppm |
| Acidity (as CH ₃ COOH) | max.0.003% |
| Assay (capillary GC) | min.99.7% |