

Code No.05-0632

# SPECIFICATION

## N-Ethyldiisopropylamine for Peptide Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Water	max.0.01%
Ethylamine (CH <sub>3</sub> CH <sub>2</sub> NH <sub>2</sub> )	max.0.0015%
N-Ethylisopropylamine [(CH <sub>3</sub> ) <sub>2</sub> CHNHC <sub>2</sub> H <sub>5</sub> ]	max.0.015%
Isopropylamine [(CH <sub>3</sub> ) <sub>2</sub> CHNH <sub>2</sub> ]	max.0.005%
Assay (capillary GC)	min.99.0%