

Code No.05-0701

SPECIFICATION

Ethanol(99.5) for Pesticide Residue And Polychlorinated Biphenyl Analysis
Concentration 5000

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
* Suitability for pesticide residue analysis (GC-FPD)	to pass test
Assay (C ₂ H ₅ OH) (GC)	min.99.5%(area/area)
Density (20°C)	0.789~0.791g/mL
Water	max.0.3%(mass/mass)
Residue after evaporation	max.2ppm(mass/mass)
Suitability for pesticide residue-PCB analysis	to pass test

JIS K 8093:2023(*:Additional test performed by Wako but not required by JIS)