

Code No.06-0539

## SPECIFICATION

### Flt3 Ligand, Human, recombinant, Animal-derived-free for Cellbiology

REQUIREMENT	SPECIFICATION
Appearance	to pass test
Bacterial endotoxins	less than 0.01 ng/ $\mu$ g
Electrophoresis test	to pass test
Mass spectrometry (MALDI-TOF)	to pass test

Source: E. coli expressed human Flt3 ligand. This product is recombinant cytokine which is produced in the process of culture and purification without using any substances from animals.