

Code No.07-0048

SPECIFICATION

D(+)-Glucuronolactone

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly brown, crystalline powder
Solubility in water	to pass test
Specific rotation [α]D ₂₀ (c=8,H ₂ O)	+18.0~+21.0°
Water	max.0.3%
Residue after ignition (as sulfate)	max.0.05%
Chloride (Cl)	max.0.001%
Sulfate (SO ₄)	max.0.005%
Copper (Cu)	max.0.001%
Lead (Pb)	max.0.001%
Iron (Fe)	max.0.001%
Assay	99.0~102.0%