

Code No.07-0050

## SPECIFICATION

## L-Glutamic Acid

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White crystalline powder
Assay (C <sub>5</sub> H <sub>9</sub> NO <sub>4</sub> )	min.99.0%(mass/mass)
Solubility in dil.hydrochloric acid	to pass test
Specific rotation [ $\alpha$ ] <sub>D20</sub>	+31.5~+32.5°
Loss on drying at 105°C	max.0.2%(mass/mass)
Residue after ignition (as sulfate)	max.0.05%(mass/mass)
Chloride (Cl)	max.0.02%(mass/mass)
Sulfate (SO <sub>4</sub> )	max.0.01%(mass/mass)
Heavy metals (as Pb)	max.0.001%(mass/mass)
Arsenic (As)	max.1ppm(mass/mass)
Iron (Fe)	max.0.001%(mass/mass)
Ammonium (NH <sub>4</sub> )	max.0.02%(mass/mass)
Foreign amino acids	to pass test

JIS K 9047:2012

\*:Additional test performed by Wako but not required by JIS