

Code No.07-0523

SPECIFICATION

Ginsenoside Rc
for the Japanese Pharmacopoeia General Tests

| REQUIREMENT | SPECIFICATION |
|--|--------------------------|
| Description (JP method) | White crystalline powder |
| Purity (JP method) | to pass test |
| * Solubility in methanol | to pass test |
| * Loss on drying over silica gel in vacuum | max.2.0% |
| * Assay (HPLC) | min.98.0% |

This product is suited for Ginsenoside Rc of general tests on The Japanese Pharmacopoeia.(*:Additional test performed by Wako)