

Code No.07-0544

SPECIFICATION

L-Glutamic Acid Standard

for Amino Acid Analysis

REQUIREMENT	SPECIFICATION
Appearance	White crystalline powder
Solubility in dil.hydrochloric acid	to pass test
Specific rotation [α]D20(c=5,dil.HCl)	+31.5~+32.5°
Loss on drying at 105°C	max.0.2%
Suitability for amino acids determination	to pass test
Assay (volumetric analysis)	min.99.0%
Assay (subtracting method)	min.99.0%