

Code No.08-0121

## SPECIFICATION

Hydroquinone  
Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White ~ nearly white, crystals
Assay	min.99.0%
Solubility in dil.acetic acid	to pass test
Melting point	171~174°C
Residue after ignition (as sulfate)	max.0.02%
Chlorine compounds (as Cl)	max.0.001%
Sulfate (SO <sub>4</sub> )	max.0.01%
Lead (Pb)	max.0.001%
Iron (Fe)	max.0.001%
Catechol (C <sub>6</sub> H <sub>6</sub> O <sub>2</sub> )	max.0.3%

JIS K 8738-1994

\*:Additional test performed by Wako but not required by JIS