

Code No.08-0715

## SPECIFICATION

### 0.5mol/l 1-Heptanesulfonic Acid Sodium Salt Solution

for Ion-pair Chromatography

| REQUIREMENT            | SPECIFICATION          |
|------------------------|------------------------|
| Appearance             | Colorless clear liquid |
| Absorbance (210nm)     | max.0.30               |
| Absorbance (220~350nm) | max.0.10               |
| pH at 25°C             | 5.0~7.0                |
| Fluorescence test      | to pass test           |
| Concentration          | 0.45~0.55mol/l         |