

Code No.08-0791

## SPECIFICATION

## Hexane for Pesticide Residue And Polychlorinated Biphenyl Tests

Concentration 5000

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
* Suitability for pesticide residue analysis (GC-FPD)	to pass test
Assay [CH <sub>3</sub> (CH <sub>2</sub> ) <sub>4</sub> CH <sub>3</sub> ](GC)	min.96.0%(area/area)
Density (20°C)	0.658~0.662g/mL
Water	max.0.01%(mass/mass)
** Residue after evaporation	max.1ppm(mass/mass)
Suitability for pesticide residue-PCB analysis	to pass test

JIS K 8825:2020(\*:Additional test performed by Wako but not required by JIS)(\*\*:The specification is originally set by Wako not JIS)