

Code No.08-1008

SPECIFICATION

Hydrogen Peroxide  
for Isolator

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Residue after evaporation	max.10ppm
Chloride (Cl)	max.1ppm
Phosphate (PO4)	max.1ppm
Sulfate (SO4)	max.1ppm
Sodium (Na)	max.50ppb
Potassium (K)	max.50ppb
Copper (Cu)	max.10ppb
Aluminium (Al)	max.500ppb
Tin (Sn)	measured value
Lead (Pb)	max.100ppb
Arsenic (As)	max.5ppb
Iron (Fe)	max.50ppb
Total nitrogen (N)	max.2ppm
Isolater test	to pass test
Assay	34.0~36.0%