

Code No.08-1045

SPECIFICATION

Hyaluronic Acid Sodium Salt, from Rooster Comb
for Biochemistry

REQUIREMENT	SPECIFICATION
Appearance	White ~ pale yellow, crystalline powder ~ powder or granules to pass test max.15.0% 15.0~20.0% 5.0~8.5 2.8~4.0%
Solubility in water	
Loss on drying at 60℃ over phosphorus(V) oxide in vacuum	
Residue after ignition (as sulfate) (after drying)	
pH of a 1g/L solution at 25℃	
Nitrogen compounds (as N)	