

SPECIFICATION

Lead(II) Nitrate

Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|--|--------------------------------------|
| * Appearance | White, crystals ~ crystalline powder |
| Assay [Pb(NO ₃) ₂] | min.99.5%(mass/mass) |
| Solubility in water | to pass test |
| Chloride (Cl) | max.0.001%(mass/mass) |
| Sodium (Na) | max.0.005%(mass/mass) |
| Potassium (K) | max.0.005%(mass/mass) |
| Copper (Cu) | max.0.001%(mass/mass) |
| Calcium (Ca) | max.0.005%(mass/mass) |
| Iron (Fe) | max.0.001%(mass/mass) |

JIS K 8563 : 2018

(*:Additional test performed by Wako but not required by JIS)