

Code No.12-0104

## SPECIFICATION

### Linseed Oil

Wako 1st Grade

| REQUIREMENT           | SPECIFICATION  |
|-----------------------|--|
| Appearance            | Slightly yellow ~ yellowish brown, clear ~ slightly muddy liquid |
| Density (20°C)        | 0.920~0.935g/ml  |
| Refractive index 20°C | 1.480~1.484  |
| Acid value            | max.1.0  |
| Saponification value  | 188~196  |
| Iodine value          | 170~204  |