

Code No.12-0437

SPECIFICATION

Liquid Paraffin

Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|--|------------------------|
| * Appearance | Colorless clear liquid |
| Density (20°C) | 0.86~0.89g/ml |
| Acidity (as H ₂ SO ₄) | max.0.002%(mass/mass) |
| Alkalinity (as NaOH) | max.0.001%(mass/mass) |
| Chloride (Cl) | max.0.001%(mass/mass) |
| Sulfate (SO ₄) | max.0.001%(mass/mass) |
| Sulfide | to pass test |
| Solid paraffin | to pass test |
| Substances darkened by sulfuric acid | to pass test |
| Substances reducing permanganate (as O) | to pass test |

JIS K 9003:2014

(*:Additional test performed by Wako but not required by JIS)