

Code No.12-0579

SPECIFICATION

Lead(II) Acetate Trihydrate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White crystals
Solubility in dil.acetic acid	to pass test
Chloride (Cl)	max.3ppm
Nitrate	to pass test (NO ₃ :max.about 0.003%)
Sodium (Na)	max.0.005%
Potassium (K)	max.0.005%
Copper (Cu)	max.5ppm
Magnesium (Mg)	max.0.005%
Calcium (Ca)	max.0.005%
Iron (Fe)	max.5ppm
Assay	min.99.0%