

Code No.12-0680

## SPECIFICATION

Locked Nucleic Acid-T Cyanoethyl Phosphoramidite (mixture of isomers)  
for Nucleic Acid Synthesis

| REQUIREMENT                       | SPECIFICATION                               |
|-----------------------------------|---------------------------------------------|
| Appearance                        | White ~ yellow, crystalline powder ~ powder |
| Solubility in acetonitrile        | to pass test                                |
| Water                             | max.0.5%                                    |
| Assay (HPLC) (mixture of isomers) | min.98.0%                                   |