

SPECIFICATION

Dichloromethane

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Description Assay (CH ₂ Cl ₂) (GC) Appearance Density (20°C) Water Residue after evaporation Acidity (as HCl) Free chlorine Iron (Fe) Substances darkened by sulfuric acid Iodine-consuming substances Stabilizer 2-Methyl-2-butene [CH ₃ CH=C(CH ₃) ₂](GC) Stabilizer Methanol (CH ₃ OH)(GC)	Colorless clear liquid min.99.5%(mass/mass) max.10 Hazen units 1.320~1.330g/mL max.0.05%(mass/mass) max.0.002%(mass/mass) max.0.01%(mass/mass) to pass test max.1ppm(mass/mass) to pass test to pass test max.0.005%(mass/mass) max.0.7%(mass/mass)

JIS K 8161:2015(*:Additional test performed by Wako but not required by JIS)