

SPECIFICATION

p-Methylstyrene

Wako 1st Grade

| REQUIREMENT | SPECIFICATION |
|-----------------------|---------------------------------------|
| Appearance | Colorless ~ pale yellow, clear liquid |
| Solubility in ethanol | to pass test |
| Density (20°C) | 0.892~0.898g/mL |
| Refractive index 20°C | 1.540~1.545 |
| Water | max.0.1% |
| o - Methylstyrene | max.3% |
| Assay (capillary GC) | min.95.0% |