

## SPECIFICATION

Multielement Standard Solution L-I  
for ICP Analysis

REQUIREMENT	SPECIFICATION
Appearance	Slightly bluish purple~pale bluish purple, clear liquid
Acidity (HNO <sub>3</sub> )	0.90~1.10mol/l
Aluminium (Al)	950~1050mg/l
Barium (Ba)	95~105mg/l
Calcium (Ca)	950~1050mg/l
Chromium (Cr)	95~105mg/l
Iron (Fe)	950~1050mg/l
Magnesium (Mg)	95~105mg/l
Lead (Pb)	95~105mg/l
Strontium (Sr)	95~105mg/l