

# SPECIFICATION

## Multication Standard Solution II

for Ion Chromatography

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Acidity (as HNO <sub>3</sub> )	0.018~0.022mol/l
Sodium ion (Na <sup>+</sup> )	95.0~105.0mg/l
Ammonium ion (NH <sub>4</sub> <sup>+</sup> )	95.0~105.0mg/l
Potassium ion (K <sup>+</sup> )	95.0~105.0mg/l