

# SPECIFICATION

## Multianion Standard Solution III for Ion Chromatography

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Bromide ion (Br <sup>-</sup> )	47.5~52.5mg/L
Chloride ion (Cl <sup>-</sup> )	9.5~10.5mg/L
Nitrite ion (NO <sub>2</sub> <sup>-</sup> )	47.5~52.5mg/L
Nitrate ion (NO <sub>3</sub> <sup>-</sup> )	47.5~52.5mg/L
Sulfate ion (SO <sub>4</sub> <sup>2-</sup> )	95~105mg/L