

Code No.14-0825

# SPECIFICATION

## Sodium Nifurstyrenate Standard for HPLC

| REQUIREMENT                             | SPECIFICATION                                                                                           |
|-----------------------------------------|---------------------------------------------------------------------------------------------------------|
| Appearance                              | Yellow ~ brown, crystalline powder ~ powder<br>or mass<br>to pass test<br><br>max.0.5%<br><br>min.98.0% |
| Solubility in sodium hydroxide solution |                                                                                                         |
| Loss on drying at 105°C                 |                                                                                                         |
| Assay (HPLC)                            |                                                                                                         |