

Code No.16-0047

# SPECIFICATION

## Liquid Paraffin

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.85~0.89g/ml
Acidity and alkalinity	to pass test
Chloride (Cl)	max.0.0015%
Sulfate (SO <sub>4</sub> )	max.0.002%
Sulfide	to pass test
Solid paraffin	to pass test
Substances darkened by sulfuric acid	to pass test
Substances reducing permanganate	to pass test