

Code No.16-0078

SPECIFICATION

Petroleum Benzine

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
Boiling range (50°C~80°C)	min.90%(v/v)
Water	max.0.05%(mass/mass)
Residue after evaporation	max.0.001%(mass/mass)
Acidity (as HCl)	max.5ppm(mass/mass)
Sulfide and reducing substances	to pass test
Fat	to pass test
Substances darkened by sulfuric acid	to pass test
* Benzene (GC)	max.1%
* Hexane (GC)	60~70%

JIS K 8594:2015

(*:Additional test performed by Wako but not required by JIS)