

Code No.16-0404

# SPECIFICATION

## Potassium Nitrate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder or mass
Solubility in water	to pass test
pH of a 50g/l solution at 25°C	5.0~8.0
Chloride (Cl)	max.0.005%
Sulfate (SO <sub>4</sub> )	max.0.01%
Copper (Cu)	max.0.001%
Sodium (Na)	max.0.02%
Lead (Pb)	max.0.002%
Iron (Fe)	max.0.001%
Assay	min.99.0%