

Code No.16-1178

SPECIFICATION

(R)-(+)-1-Phenylethylamine

Wako Special Grade

| REQUIREMENT | SPECIFICATION |
|---|---|
| Appearance | Colorless ~ slightly yellow, clear liquid |
| Density (20°C) | 0.948~0.956g/ml |
| Refractive index 20°C | 1.524~1.529 |
| Specific rotation [α]D20(neat) | +37~+41° |
| Water | max.1.0% |
| Assay (capillary GC) | min.98.0% |