

Code No.16-1199

# SPECIFICATION

Propiophenone

Wako Special Grade

| REQUIREMENT           | SPECIFICATION                         |
|-----------------------|---------------------------------------|
| Appearance            | Colorless ~ pale yellow, clear liquid |
| Density (20°C)        | 1.007~1.013g/mL                       |
| Freezing point        | 17~20°C                               |
| Refractive index 20°C | 1.524~1.527                           |
| Water                 | max.0.2%                              |
| Assay (capillary GC)  | min.99.0%                             |