

Code No.16-1701

SPECIFICATION

Pyridine

Infinity Pure

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Solubility in water	to pass test
Water	max.0.1%
Residue after evaporation	max.20ppm
Chloride (Cl)	max.10ppm
Sulfate (SO ₄)	max.10ppm
Copper (Cu)	max.5ppb
Lead (Pb)	max.5ppb
Iron (Fe)	max.10ppb
Ammonium (NH ₄)	max.20ppm
Peroxides (as H ₂ O ₂)	max.10ppm
Substances reducing permanganate (as O)	max.20ppm
Tollens reactive substances (as CH ₃ CHO)	max.10ppm
Assay (capillary GC)	min.99.5%