

Code No.16-1791

SPECIFICATION

Polyvinyl Alcohol

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White ~ slightly yellow,shot ~ powder
Solubility in water	to pass test
Loss on drying at 105°C	max.10.0%(mass/mass)
Residue after ignition (as sulfate)	max.1.0%(mass/mass)
Viscosity (4%, 20°C)	35~55mPa·s
Suitability for absorptiometric determination	to pass test
Saponification degree	78~82mol%

JIS K 9550 : 2020(*:Additional test performed by Wako but not required by JIS)