

Code No.16-1792

## SPECIFICATION

## Phenylfluorone

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Yellowish red ~ brown, crystalline powder ~ powder to pass test
Solubility in dil. hydrochloric acid-methanol	
Absorbance (7.5mg/L)	min.0.70
Loss on drying at 105°C	max.5.0%(mass/mass)
Residue after ignition (as sulfate)	max.0.5%(mass/mass)
Suitability for tin determination	to pass test

JIS K 9547 : 2018

(\*:Additional test performed by Wako but not required by JIS)