

Code No.16-2279

SPECIFICATION

Polysorbate 20

Japanese Pharmaceutical Excipients

REQUIREMENT	SPECIFICATION
* Description	Slightly yellow ~ yellow, liquid
Identification	to pass test
Viscosity	350~550mm2/s
pH	4.0~7.0
Sp. Gr. d 2020	1.090~1.130
Acid value	not more than 4.0
Saponification value	43~55
Purity	-
(1) Heavy metals	not more than 20ppm
(2) Arsenic	not more than 2ppm
Water	not more than 3.0%
Residue on ignition	not more than 0.25%

Japanese Pharmaceutical Excipients(\*:Reference item)