

Code No.16-2367

SPECIFICATION

Polysorbate 80 Standard

for Food Additives Test

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ pale yellow, clear ~ slightly muddy liquid
Kinematic viscosity (25°C)	345~445mm ² /s
Sp. Gr.(20/20°C)	1.065~1.095
Refractive index 20°C	under investigation
Water	max.3.0%
Acid value	max.2
Saponification value	45~55
Iodine value	19~24
Hydroxyl value	under investigation
Assay (capillary GC)	min.98.0%