

SPECIFICATION

Palladium(II) Nitrate

for Organic Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Brown ~ blackish brown, powder or mass or granules
Solubility in dil.hydrochloric acid	to pass test
Sodium (Na)	max.0.05%
Magnesium (Mg)	max.0.005%
Aluminium (Al)	max.0.005%
Potassium (K)	max.0.5%
Calcium (Ca)	max.0.05%
Manganese (Mn)	max.0.005%
Iron (Fe)	max.0.005%
Nickel (Ni)	max.0.005%
Cobalt (Co)	max.0.005%
Copper (Cu)	max.0.005%
Zinc (Zn)	max.0.01%
Ruthenium (Ru)	max.0.01%
Rhodium (Rh)	max.0.05%
Iridium (Ir)	max.0.05%
Barium (Ba)	max.0.005%
Lead (Pb)	max.0.05%
Platinum (Pt)	max.0.01%
Assay	min.98.0%