

SPECIFICATION

2-Propanol
for Medicine Test

REQUIREMENT	SPECIFICATION
* Appearance Assay [(CH ₃) ₂ CHOH] (GC) Solubility in water Density (20°C) Refractive index 20°C Water Residue after evaporation Acidity (as C ₂ H ₅ COOH) Carbonyl compounds (as CO) Methanol (CH ₃ OH) (GC) Iron (Fe) Substances darkened by sulfuric acid Sp. Gr.(20/20°C) (Ph.Eur.)	Colorless clear liquid min.99.7%(mass/mass) to pass test 0.784~0.786g/ml 1.376~1.379 max.0.1%(mass/mass) max.0.001%(mass/mass) max.0.002%(mass/mass) max.0.005%(mass/mass) max.0.1%(mass/mass) max.0.1 ppm(mass/mass) to pass test about 0.785

JIS K 8839-2007

*:Additional test performed by Wako but not required by JIS