

Code No.16-2484

SPECIFICATION

1-Propanol, Super Dehydrated for Organic Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Solubility in water	to pass test
Density (20°C)	0.801~0.806g/ml
Refractive index 20°C	1.383~1.387
Water	max.0.001%
Residue after evaporation	max.0.001%
Acidity (as C ₂ H ₅ COOH)	max.0.002%
Substances darkened by sulfuric acid	to pass test
Assay (capillary GC)	min.99.5%