

Code No.16-2485

## SPECIFICATION

2-Propanol, Super Dehydrated  
for Organic Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Solubility in water	to pass test
Density (20°C)	0.784~0.787g/ml
Refractive index 20°C	1.376~1.379
Water	max.0.001%
Residue after evaporation	max.0.001%
Acidity (as C <sub>2</sub> H <sub>5</sub> COOH)	max.0.002%
Substances darkened by sulfuric acid	to pass test
Assay (capillary GC)	min.99.7%