

SPECIFICATION

Paeonol

for the Japanese Pharmacopoeia Crude Drugs Test (for Assay and TLC)

| REQUIREMENT | SPECIFICATION |
|---|--------------------------------------|
| Description (JP method) | White crystals or crystalline powder |
| Purity (JP method) | - |
| Related substances (JP method) | to pass test |
| Unity of peak (JP method) | to pass test |
| * Assay (qNMR) (JP method) | min.98.0% |
| ** Loss on drying over calcium chloride | max.1.0% |

This product is suited for Paeonol for assay and thin-layer chromatography of general tests on The Japanese Pharmacopoeia.
(*:The specification is originally set by Wako)(**):Additional test performed by Wako)