

Code No.16-2494

SPECIFICATION

Polyethylene Glycol 4,000,000

Wako 1st Grade

| REQUIREMENT | SPECIFICATION |
|-------------------------------------|---|
| Appearance | White ~ slightly yellow, crystalline powder ~ |
| Viscosity (0.5%, 25°C) | powder or granules under investigation |
| Solubility in water | to pass test |
| Loss on drying at 105°C | max.2.0% |
| Residue after ignition (as sulfate) | max.2.0% |
| pH of a 5g/l solution at 25°C | 6.0~8.0 |
| Acidity (as CH ₃ COOH) | max.0.2% |