

Code No.18-0266

SPECIFICATION

Rosmarinic Acid

for the Japanese Pharmacopoeia Crude Drugs Test (for Assay and Thin-Layer Chroma

REQUIREMENT	SPECIFICATION
Description (JP method)	White to pale yellow, crystals or crystalline powder
Identification (JP method)	to pass test
Purity (JP method)	-
Related substances (JP method)	to pass test
Unity of peak (JP method)	to pass test
* Assay (qNMR) (JP method)	min.95.0%
** Solubility in ethanol (after drying)	to pass test
** Loss on drying over silica gel in vacuum	max.1.0%

This product is suited for Rosmarinic Acid for assay and thin-layer chromatography of general tests on The Japanese Pharmacopoeia.(*:The specification is originally set by Wako)(**:.Additional test performed by Wako)