

SPECIFICATION

Silver(I) Sulfate

Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|--|---|
| * Appearance Assay (Ag ₂ SO ₄) Solubility in dil.nitric acid Nitrogen compounds (as N) Substances not precipitated by HCl Copper (Cu) Iron (Fe) | White crystalline powder min.99.5%(mass/mass) to pass test max.0.001%(mass/mass) max.0.03%(mass/mass) max.0.004%(mass/mass) max.0.004%(mass/mass) |

JIS K 8965 : 2019(*:Additional test performed by Wako but not required by JIS)