

SPECIFICATION

Silver(I) Sulfate
Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|--|--------------------------|
| * Appearance | White crystalline powder |
| Assay (Ag ₂ SO ₄) | min.99.5%(mass/mass) |
| Solubility in dil.nitric acid | to pass test |
| Nitrogen compounds (as N) | max.0.001%(mass/mass) |
| Substances not precipitated by HCl | max.0.03%(mass/mass) |
| Copper (Cu) | max.0.004%(mass/mass) |
| Iron (Fe) | max.0.004%(mass/mass) |

JIS K 8965 : 2019(*:Additional test performed by Wako but not required by JIS)