

Code No.19-0137

SPECIFICATION

Sodium Hydrogensulfite

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White ~ nearly white, crystalline powder ~
Assay (as SO ₂)	64.0~67.4%(mass/mass)
Solubility in water	to pass test
Chloride (Cl)	max.0.01%(mass/mass)
Potassium (K)	max.0.01%(mass/mass)
Copper (Cu)	max.0.001%(mass/mass)
Lead (Pb)	max.0.001%(mass/mass)
Arsenic (As)	max.5ppm(mass/mass)
Iron (Fe)	max.0.002%(mass/mass)

JIS K 8059:2018(*:Additional test performed by Wako but not required by JIS)