

Code No.19-0141

SPECIFICATION

Sodium Tetraborate Decahydrate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White, crystals ~ crystalline powder
Assay ($\text{Na}_2\text{B}_4\text{O}_7 \cdot 10\text{H}_2\text{O}$)	99.5~101.0%(mass/mass)
Solubility in water	to pass test
Chloride (Cl)	max.5ppm(mass/mass)
Nitrate	to pass test
Phosphate (PO_4)	max.0.001%(mass/mass)
Sulfate (SO_4)	max.0.003%(mass/mass)
Copper (Cu)	max.5ppm(mass/mass)
Calcium (Ca)	max.0.001%(mass/mass)
Lead (Pb)	max.5ppm(mass/mass)
Arsenic (As)	max.2ppm(mass/mass)
Iron (Fe)	max.3ppm(mass/mass)

JIS K 8866:2008(*:Additional test performed by Wako but not required by JIS)