

## SPECIFICATION

## Sodium Formate

## Guaranteed Reagent

| REQUIREMENT  | SPECIFICATION  |
|--|--|
| * Appearance<br>Assay (HCOONa)<br>Solubility in water<br>pH of a 50g/l solution at 25°C<br>Chloride (Cl)<br>Phosphate (PO <sub>4</sub> )<br>Sulfate (SO <sub>4</sub> )<br>Magnesium (Mg)<br>Calcium (Ca)<br>Lead (Pb)<br>Iron (Fe) | White, crystals ~ crystalline powder<br>min.98.0%(mass/mass)<br>to pass test<br>7.0~8.5<br>max.0.001%(mass/mass)<br>max.0.001%(mass/mass)<br>max.0.001%(mass/mass)<br>max.0.001%(mass/mass)<br>max.0.001%(mass/mass)<br>max.5ppm(mass/mass)<br>max.5ppm(mass/mass) |

JIS K 8267-2014

(\*:Additional test performed by Wako but not required by JIS)