

SPECIFICATION

Trisodium Phosphate Dodecahydrate

Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|---|--|
| * Appearance Assay (Na ₃ PO ₄ ·12H ₂ O) Solubility in water Chloride (Cl) Nitrate Sulfate (SO ₄) Calcium (Ca) Lead (Pb) Arsenic (As) Iron (Fe) Sodium hydroxide (NaOH) | White, crystals ~ crystalline powder min.99.0%(mass/mass) to pass test max.0.002%(mass/mass) to pass test max.0.005%(mass/mass) max.0.02%(mass/mass) max.0.001%(mass/mass) max.1ppm(mass/mass) max.0.001%(mass/mass) max.1.5%(mass/mass) |

JIS K 9012:2024(*:Additional test performed by Wako but not required by JIS)