

Code No.19-0300

## SPECIFICATION

## Potassium Sodium (+)-Tartrate Tetrahydrate

Wako 1st Grade

| REQUIREMENT                            | SPECIFICATION  |
|--|----------------|
| Appearance                             | White crystals |
| Solubility in water                    | to pass test   |
| pH of a 50g/l solution at 25°C         | 6.5~8.5        |
| Chloride (Cl)                          | max.0.003%     |
| Sulfur compounds (as SO <sub>4</sub> ) | max.0.03%      |
| Copper (Cu)                            | max.4ppm       |
| Lead (Pb)                              | max.4ppm       |
| Iron (Fe)                              | max.0.002%     |
| Ammonium (NH <sub>4</sub> )            | max.0.005%     |
| Assay                                  | min.99.0%      |